

IHE Integration Statement

Vendor	Product Name	Version	Date
Aspyra, Inc.	AccessNET	6.2	8/30/2006

This product is intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Scheduled Workflow	<ul style="list-style-type: none"> ^ Image Manager ^ Image Creator ^ Image Display 	None
Patient Information Reconciliation	<ul style="list-style-type: none"> ^ Image Manager 	None
Consistent Presentation of Images	<ul style="list-style-type: none"> ^ Image Manager ^ Image Creator ^ Image Display ^ Print Composer 	None
Presentation of Grouped Procedures	<ul style="list-style-type: none"> ^ Image Manager ^ Image Display 	None
Access to Radiology Information	<ul style="list-style-type: none"> ^ Image Manager ^ Image Display 	None
Key Image Notes	<ul style="list-style-type: none"> ^ Image Manager ^ Image Creator ^ Image Display 	None
Simple Image and Numeric Reports	<ul style="list-style-type: none"> ^ Report Repository ^ Image Creator ^ Image Display 	None
Portable Data for Imaging	<ul style="list-style-type: none"> ^ Portable Media Creator ^ Portable Media Importer ^ Image Display ^ Report Reader ^ Print Composer 	None

Aspyra Internet Address: www.aspyra.com

Links to Standards Conformance Statements for the Implementation

HL7 www.aspyra.com/solutions/integration

DICOM www.aspyra.com/solutions/integration/dicom_accessnet.pdf

Links to general information on IHE

In North America:
www.ihe.net

In Europe:
www.ihe-europe.org

In Japan:
www.iira-net.or.jp/ihe-j

